



Luminaire Based Accessoires

LRM8140/00 ActiLume G2 RF Mov Sensor

Radio frequency (RF) radar motion detection sensor which is to be connected to an ActiLume Wireless DALI gen2 controller.

Product data

Product Data		Numerator – packs per outer box	
Full product code	871869656919100	Numerator – packs per outer box	12
Order product name	LRM8140/00 ActiLume G2 RF Mov Sensor	SAP material	913700389103
EAN/UPC – product	8718696569191	SAP net weight (piece)	33.000 g
Order code	56919100		
SAP numerator – quantity per pack	1		